



780 BOSTON ROAD  
BILLERICA, MA 01821  
T 978.663.0212 | F 978.663.0214  
[BELMONTMEDTECH.COM](http://BELMONTMEDTECH.COM)

Subject: Material Components

March 18, 2019

Dear Customer,

Please be informed that aluminum is not used in the manufacturing of the Belmont disposables of Rapid Infuser, buddy™ and Hyperthermia pump™ product lines as listed below:

- Part number 903-00006, 3- Spike Disposable Set (12 per box)
- Part number 903-00106, 3- Spike Disposable Set (6 per box)
- Part number 903-00032, 4x4 Disposable Set ((4) LVR + (4) Heater Exchangers w/Pt. Line)
- Part number 903-00018, 3 Liter Reservoir (6 per box)
- Part number 903-00004, Dual Patient Line (18 per box)
- Part number 903-00022, Patient Line Extension (24 per box)
- Part number 905-00010, buddy™ Disposable Set (24 per case)
- Part number 902-00045, HP Procedure Kit (1 kit)
- Part number 902-00037, HP Straight Set (1 kit)
- Part number 902-00038, HP Bifurcated Set (1 kit)

Additionally, these consumable devices have been tested for biocompatibility in alignment with ISO10993-1 Biological Evaluation of Medical Devices, and the fluid path has been tested and shown to be fully biocompatible and have shown no mutagenicity nor cellular death.

For any additional technical questions, please contact our Technical Service Team at [techservice@belmontmedtech.com](mailto:techservice@belmontmedtech.com) or call +978.663.0212, extension 2.

For Sales & Customer Service, please contact us at [orders@belmontmedtech.com](mailto:orders@belmontmedtech.com) or call +978.663.0212, extension 1.

Thank you and Best Regards,

Belmont Medical Technologies