




Arrow

Paediatric CVC Sets

Nitinol guidewire and safety features

Teleflex's Arrow Paediatric CVC Sets

Although generally safe and routine in ICU, central venous catheter insertion procedures can be challenging and involve complications, especially where a child is concerned.

Teleflex's Arrow has always strived to provide the optimum quality products to ensure ease-of-use and the utmost in patient safety and comfort. That is why arrow central venous catheters are the first choice for venous access.

Arrow's paediatric CVCs are, of course, no exception. The flexible polyurethane they are made of is durable but pliant and soft enough to protect the vessel.

Arrow's paediatric CVCs come in multiple lumen and length configurations to meet the specific therapy needs of children. Available configurations range from 1 to 3 lumens and come in several lengths and French sizes. You also have the choice of uncoated or coated catheters for antimicrobial protection.

Set Components

- Polyurethane catheter (PUR) with Blue FlexTip
- SWG Nitinol guidewire and arrow Advancer
- Introducer catheter over needle
- 5 ml Luer slip syringe
- Introducer needles
- 2 dilators
- Catheter clamp and fastener
- Safety scalpel
- SharpsAway II Locking Disposal Cup

Teleflex's Arrow Paediatric CVC Sets with Nitinol Guidewire, Introducer Needle and Additional Safety Features


In order to meet customer needs, we have reviewed the content of our paediatric CVC sets and now offer the following:




Nitinol guidewire: the Nitinol guidewire has a welded connection at the tip, a high tensile strength and is highly flexible. Additionally the guidewire has cm-markings along the shaft and a small J tip radius.



New introducer needle: the introducer needle provided is siliconised and has a sharper bevel, this introducer is used in our other products. These features can provide more comfort for your patient and make it easier for the clinical staff to attend to the task at hand.



SharpsAway II Locking Disposal Cup and safety scalpel: the avoidance of needle stick injuries for your staff not only complies with key regulatory guidelines, it's being responsible by minimising the risk to your patient and the clinical staff delivering life-saving care. Our safety scalpel is easy-to-control to avoid sharps injuries, while our SharpsAway II Locking Disposal Cup allows safe disposal of used sharps to minimise the risk of exposure.



Paediatric Nitinol Guidewire Kits – with cm-marking, Latex-free

REF.	CATHETER	CROSS SECTION	LUMEN	DIAMETER	TIP	MATERIAL	QTY
EU-04522-N	22 Ga×4" (10 cm)	n/a	1	.021" (0.53 mm) × 13¾" (35 cm)	Straight/Straight Nitinol	PUR	5
EU-04150-N	20 Ga×5" (12 cm)	n/a	1	.025" (0.64 mm) × 13¾" (35 cm)	J/straight Nitinol	PUR	5
EU-04218-N	18 Ga×8" (20 cm)	n/a	1	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-12402-N	4 Fr×2" (5 cm)	C-22/22	2	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-15402-N	4 Fr×3⅛" (8 cm)	C-22/22	2	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-25402-N	4 Fr×3⅛" (8 cm)	C-22/22	2	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR/AGB	5
EU-14402-N	4 Fr×5⅛" (13 cm)	C-22/22	2	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-24402-N	4 Fr×5⅛" (13 cm)	C-22/22	2	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR/AGB	5
EU-16402-N	4 Fr×12" (30 cm)	C-22/22	2	.018" (0.46 mm) × 26¾" (68 cm)	J/straight Nitinol	PUR	5
EU-14502-N	5 Fr×5⅛" (13 cm)	E-18/20	2	.021" (0.53 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-22502-N	5 Fr×2" (5 cm)	E-18/20	2	.021" (0.53 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR/AGB	5
EU-15403-N	4 Fr×3" (8 cm)	H-20/22/22	3	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-14403-N	4 Fr×5⅛" (13 cm)	H-20/23/23	3	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-15553-N	5,5 Fr×3" (8 cm)	H-20/22/22	3	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-25553-N	5,5 Fr×3⅛" (8 cm)	H-20/22/22	3	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR/AGB	5
EU-16553-N	5,5 Fr×5⅛" (13 cm)	H-20/22/22	3	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR	5
EU-26553-N	5,5 Fr×5⅛" (13 cm)	H-20/22/22	3	.018" (0.46 mm) × 17¾" (45 cm)	J/straight Nitinol	PUR/AGB	5
EU-14553-N	5,5 Fr×12" (30 cm)	H-20/22/22	3	.018" (0.46 mm) × 26¾" (68 cm)	J/straight Nitinol	PUR	5

AGB = ARROWg^{ard} Blue

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose-driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow, Deknatel, Hudson RCI, LMA, Pilling, QuickClot, Rüsch, UroLift, and Weck – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia 1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +1 (0)800 387 9699

China (Shanghai) +86 (0)21 6060 7000

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 691 3004

Greece +30 210 67 77 717

India +91 (0)80 4093 4790

Italy +39 0362 58911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand 0800 601 100

Poland +48 22 4624032

Portugal +351 225 491 051

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 91 198 84 31

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com

Teleflex, the Teleflex logo, Arrow, SharpsAway and Advancer are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 05/2021. Subject to technical changes without further notice. © 2021 Teleflex Incorporated. All rights reserved.